

DigiTUBEs

Disposable digestion tubes

DigiTUBEs are disposable digestion tubes suitable for use with commercially available graphite digestion block systems.

DigiTUBEs are manufactured with ICP-MS Grade Polypropylene and specially designed with a flat bottom to provide more direct contact with the base of digestion block wells, resulting in a faster heat-up time and more efficient sample temperature response times. Every case of **DigiTUBEs** comes complete with a lot specific Calibration Certificate that guarantees Class A tolerances for every **DigiTUBE**.



Products meet the following tolerance requirements of ASTM Standard E1272:

- ▶ ± 0.10 ml for the 15 ml volume line on 15 ml tubes (Class A)
- ▶ ± 0.25 ml for all volume lines on 50 ml tubes (Class A)
- ▶ ± 0.50 ml for 50 ml and 100 ml volume lines on 100 ml tubes (Class A)

For clean applications, **DigiTUBEs** are also available in shrink-wrapped bundles to protect against dust infiltration during storage and handling.

Key Features

- ▶ Ultra-low leachable metal content to eliminate contamination risk
- ▶ Optional RackLock feature enables one-handed closure of screw caps
- ▶ Class A graduations for easy and accurate dilution
- ▶ Disposable, eliminating the need to clean vessel
- ▶ Leak-free screw cap available in variety of colors

DigiTUBEs can replace 5 vessels in the lab - digestion beaker, graduated cylinder, volumetric flask, autosampler tube, and storage container. Save time and lower costs with the use of **DigiTUBE**s.

Bulk DigiTUBEs

Description	Quantity	Part Number
15 ml with RackLock, blue caps	1620	010-515-261
15 ml with RackLock, no caps	1620	010-515-262
15 ml without RackLock, blue caps	1620	010-515-263
15 ml without RackLock, no caps	1620	010-515-264
15 ml with RackLock, blue caps	540	010-515-361
15 ml without RackLock, no caps	540	010-515-363
50 ml with RackLock, blue caps	750	010-500-261
50 ml with RackLock, no caps	750	010-500-262
50 ml without RackLock, blue caps	750	010-500-263
50 ml without RackLock, no caps	750	010-500-264
100 ml with RackLock, blue caps	300	010-501-261
100 ml with RackLock, no caps	300	010-501-262
100 ml without RackLock, blue caps	300	010-501-263
100 ml without RackLock, no caps	300	010-501-264

Bulk Caps

	Description	Quantity	Part Number
Screw Caps	Screw Cap, Blue 15 ml DigiTUBE	540	010-515-060
	Screw Cap, Clear, 15 ml DigiTUBE	540	010-515-172
Color Coded 50 ml DigiTUBE s	Blue Caps for 50 ml DigiTUBE s	250	010-500-060
	Yellow Caps for 50 ml DigiTUBE s	250	010-500-130
	Orange Caps for 50 ml DigiTUBE s	250	010-500-140
	Red Caps for 50 ml DigiTUBE s	250	010-500-150
	Clear Caps for 50 ml DigiTUBE s	250	010-500-172
Color Coded 100 ml DigiTUBE s	Blue Caps for 100 ml DigiTUBE s	100	010-501-060

HandyPAKs

HandyPAKs include multiple racks of **DigiTUBE**s (50 ml or 100 ml), shrink wrapped to eliminate contamination during storage or transport. Caps are provided in a separate bag. The HandyPAK Plus option contains pre-capped **DigiTUBE**s.

	Description	Quantity	Part Number
HandyPAK	50 ml DigiTUBE with RackLock, blue caps	30 x 25 tubes	010-500-362
	50 ml DigiTUBE without RackLock, blue caps	30 x 25 tubes	010-500-364
	100 ml, DigiTUBE with RackLock, blue caps	9 x 34 tubes	010-501-362
	100 ml, DigiTUBE without RackLock, blue caps	9 X 34 tubes	010-501-364
HandyPAK Plus	50 ml DigiTUBE with RackLock, blue caps	20 x 24 tubes	010-500-461
	50 ml DigiTUBE without RackLock, blue caps	20 x 24 tubes	010-500-463
	100 ml, DigiTUBE with RackLock, blue caps	24 x 12 tubes	010-501-462
	100 ml, DigiTUBE without RackLock, blue caps	24 x 12 tubes	010-501-464

Watch Glasses & Storage Racks

Description	Quantity	Part Number
Disposable Polypropylene Watch Glasses, 50 ml	1000	010-500-081
Disposable Polypropylene Watch Glasses, 100 ml	500	010-501-081
Plastic Storage Racks for 50 ml DigiTUBE	5	010-500-026
Foam Storage Racks for 100 ml DigiTUBE	5	010-501-028

Digestion Rack for 15 ml, 50 ml and 100 ml **DigiTUBE**s available.